



IMPLANTE DE MENTÓN ANATÓMICO




Referencia	Tamaño	A	B	C	D	E	F
EAC-S	Pequeño	5.5	1.0	0.45	9.0	3.7	0.2
EAC-M	Mediano	5.6	1.2	0.50	9.5	4.2	0.4
EAC-L	Grande	6.8	1.4	0.60	11.0	4.8	0.4
EAC-XL	Extragrande	6.0	1.4	0.75	11.0	5.1	0.5

IMPLANTE MALAR




Referencia	Tamaño	A	B	C
BSM-S	Pequeño	4.0	1.6	0.45
BSM-M	Mediano	4.4	1.8	0.50
BSM-L	Grande	4.9	2.0	0.57
BSM-XL	Extragrande	5.3	2.2	0.61

IMPLANTE NASAL





Referencia	Tipo	A	B	C	D	E	F	G
DCI-03S	Suave	5.8	3.0	1.1	0.3	0.9	0.8	0.6

IMPLANTE MANDIBULAR POSTERIOR



Referencia	A	B	C	D	E	F	G
TMA-6L	4.5	0.6	4.5	3.5	3.0	0.36	0.4
TMA-6R	4.5	0.6	4.5	3.5	3.0	0.36	0.4
TMA-10L	4.5	1.0	4.5	3.5	3.0	0.36	0.4
TMA-10R	4.5	1.0	4.5	3.5	3.0	0.36	0.4

BLOQUE DE SILICONA (Consistencia blanda y media)



Referencia	A	B	C
BLK-10	12.0	7.0	1.5
BLK-20	12.0	7.0	1.5